

220: CASO AGAR (Merck 105458)

This recipe contains strain-specific modifications for *Paenibacillus phyllosphaerae* DSM 17399 *

Final pH: 7.3

Final volume: 1000 ml

Casein peptone	15.00	g
Soy peptone	5.00	g
NaCl	5.00	g
Agar	15.00	g
Xylan	7.00	g/l
Yeast	4.00	g/l
MnSO ₄	10.00	mg/l
Distilled water	1000.00	ml

Adjust pH to 7.3. Medium is identical with Tryptone Soya Agar (Oxoid CM 131).

* Plus 0.7% xylan plus 0.4% yeast plus 10 µg/ml MnSO₄