

220: CASO AGAR (Merck 105458)

This recipe contains strain-specific modifications for *Amphibacillus sediminis* DSM 21624 *

Final pH: * 8.5 - 9.0

Final volume: 1000 ml

Casein peptone	15.00	g
Soy peptone	5.00	g
NaCl	5.00	g
Agar	15.00	g
Distilled water	1000.00	ml

Adjust pH to 7.3. Medium is identical with Tryptone Soya Agar (Oxoid CM 131).

* Adjust pH of final medium to 8.5 - 9.0.