

220: CASO AGAR (Merck 105458)

This recipe contains strain-specific modifications for *Paenibacillus campinasensis* DSM 21989 *

Final pH: * 10.0

Final volume: 1000 ml

Casein peptone	15.00	g
Soy peptone	5.00	g
NaCl	5.00	g
Agar	15.00	g
Distilled water	1000.00	ml

Adjust pH to 7.3. Medium is identical with Tryptone Soya Agar (Oxoid CM 131).

* Plus Na-sesquicarbonate solution (see medium 31); pH 10