

220: CASO AGAR (Merck 105458)

This recipe contains strain-specific modifications for *Paenibacillus residui* DSM 22072 *

Final pH: * 7.2

Final volume: 1000 ml

Casein peptone	15.00	g
Soy peptone	5.00	g
NaCl	5.00	g
Agar	15.00	g
MnSO ₄	10.00	µg/ml
Distilled water	1000.00	ml

Adjust pH to 7.3. Medium is identical with Tryptone Soya Agar (Oxoid CM 131).

* Plus 10 µg/ml MnSO₄; pH 7.2