

220: CASO AGAR (Merck 105458)

This recipe contains strain-specific modifications for *Geomicrobium sediminis* DSM 23309 *

Final pH: 7.3

Final volume: 1000 ml

Casein peptone	15.00	g
Soy peptone	5.00	g
NaCl	5.00	g
Agar	15.00	g
Sodium sesquicarbonate	10.00	%
Distilled water	1000.00	ml

Adjust pH to 7.3. Medium is identical with Tryptone Soya Agar (Oxoid CM 131).

* plus 10% sodium sesquicarbonate