

220: CASO AGAR (Merck 105458)

This recipe contains strain-specific modifications for *Cytobacillus purgationiresistens* DSM 23494 *

Final pH: * 7.0

Final volume: 1000 ml

Casein peptone	15.00	g
Soy peptone	5.00	g
NaCl	5.00	g
Agar	15.00	g
MnSO ₄	10.00	mg/L
Distilled water	1000.00	ml

Adjust pH to 7.3. Medium is identical with Tryptone Soya Agar (Oxoid CM 131).

* Plus 10 mg/L MnSO₄; pH 7-10