

220: CASO AGAR (Merck 105458)

This recipe contains strain-specific modifications for *Methylobacillus* sp. DSM 5746 *

Final pH: 7.3

Final volume: 1000 ml

Casein peptone	15.00	g
Soy peptone	5.00	g
NaCl	5.00	g
Agar	15.00	g
Methanol	10.00	ml/l
Distilled water	1000.00	ml

Adjust pH to 7.3. Medium is identical with Tryptone Soya Agar (Oxoid CM 131).

* Supplement medium with 10 ml/l Methanol.