

220: CASO AGAR (Merck 105458)

This recipe contains strain-specific modifications for *Pseudaminobacter salicylatoxidans* DSM 7005 *

Final pH: 7.3

Final volume: 1000 ml

| | | |
|-----------------|---------|----|
| Casein peptone | 15.00 | g |
| Soy peptone | 5.00 | g |
| NaCl | 5.00 | g |
| Agar | 15.00 | g |
| Distilled water | 1000.00 | ml |

Adjust pH to 7.3. Medium is identical with Tryptone Soya Agar (Oxoid CM 131).

* Rehydrate and grow lyophilized cells from the ampoule in 5 ml liquid broth 220. Subsequent subculturing may be carried out in liquid medium or on agar plates and in medium 465 plus gentisinic acid