

## 220: CASO AGAR (Merck 105458)

This recipe contains strain-specific modifications for *Pseudomonas fluorescens* DSM 7355 \*

Final pH: 7.3

Final volume: 1000 ml

Casein peptone	15.00	g
Soy peptone	5.00	g
NaCl	5.00	g
Agar	15.00	g
Distilled water	1000.00	ml

Adjust pH to 7.3. Medium is identical with Tryptone Soya Agar (Oxoid CM 131).

\* Rehydrate and grow lyophilized cells from the ampoule in 5 ml liquid broth 220. Subsequent subculturing may be carried out in liquid medium or on agar plates