Microorganisms



360: YPM MEDIUM

This recipe contains strain-specific modifications for Frateuria aurantia DSM 6220 *

Final volume: 1000 ml

Yeast extract	5.00	g
Peptone	3.00	g
Mannitol	25.00	g
Distilled water	1000.00	ml
Agar, for solid medium	12.00	g/l

1. For solid media add agar at a concentration of 12.0 g/l

2. pH not adjusted.

* Rehydrate and grow cells from freeze dried ampoule in 5 ml liquid medium 360. Subsequent subcultivation may be carried out in liquid medium or on agar plates