

535: TRYPTICASE SOY BROTH AGAR

This recipe contains strain-specific modifications for *Vogesella mureinivorans* DSM 21247 *

Final pH: 7.3

Final volume: 1000 ml

Trypticase soy broth (BBL 11768, Oxoid CM129 or Merck 30460)	30.00	g
Agar	20.00	g
Distilled water	1000.00	ml

1. pH 7.3
2. Autoclave at 121°C for 15 min.

* Rehydrate and grow lyophilized cells from the ampoule in 5 ml liquid broth 545.
Subsequent subculturing may be carried out in liquid medium 545 or on agar plates