

535: TRYPTICASE SOY BROTH AGAR

This recipe contains strain-specific modifications for *Pseudomonas fluorescens* DSM 4358 *

Final pH: 7.3

Final volume: 1000 ml

Trypticase soy broth (BBL 11768, Oxoid CM129 or Merck 30460)	30.00	g
Agar	20.00	g
Distilled water	1000.00	ml

1. pH 7.3
2. Autoclave at 121°C for 15 min.

* For fluorescence; Rehydrate and grow lyophilized cells from the ampoule in 5 ml liquid broth 1. Subsequent subculturing may be carried out in liquid medium or on agar plates