

535: TRYPTICASE SOY BROTH AGAR

This recipe contains strain-specific modifications for *Pseudomonas mosselii* DSM 17497 *

Trypticase soy broth (BBL 11768, Oxoid CM129 or Merck 30469)	30.00	g
Agar	20.00	g
Distilled water	1000.00	ml

1. pH 7.3
2. Autoclave at 121°C for 15 min.

* Rehydrate and grow lyophilized cells from the ampoule in 5 ml liquid broth. Subsequent subculturing may be carried out in liquid medium or on agar plates