

## 535: TRYPTICASE SOY BROTH AGAR

This recipe contains strain-specific modifications for *Photorhabdus temperata* DSM 12190 \*

Final pH: 7.3

Final volume: 1000 ml

Trypticase soy broth (BBL 11768, Oxoid CM129 or Merck <del>30460</del> )	30.00	g
Agar	20.00	g
Distilled water	1000.00	ml

1. pH 7.3
2. Autoclave at 121°C for 15 min.

\* Resuspend and grow freeze dried cells in liquid medium. Subsequent subculturing may be carried out on plates or in liquid medium.