

## 53: CORYNEBACTERIUM AGAR

This recipe contains strain-specific modifications for *Marinococcus halotolerans* DSM 16375 \*

Final pH: 7.2 - 7.4

Final volume: 1000 ml

Casein peptone (tryptic digest)	10.00	g
Yeast extract	5.00	g
Glucose	5.00	g
NaCl	5.00	g
Agar	15.00	g
MgCl <sub>2</sub>	100.00	g/l
Distilled water	1000.00	ml

Adjust pH to 7.2 - 7.4.

\* Plus 10% MgCl<sub>2</sub>