

## 636: ISP MEDIUM 2 WITH 5% NaCl

This recipe contains strain-specific modifications for *Paracoccus saliphilus* DSM 18447 \*

Final pH: 7.3 - 7.5

Final volume: 1000 ml

Yeast extract (Difco)	4.00	g
Malt extract (Difco)	10.00	g
Dextrose (Difco)	4.00	g
Agar	20.00	g
NaCl (5 %)	20.00	g
Distilled water	1000.00	ml

pH 7.3 - 7.5

\* Reduce the amount of NaCl to 2%