

## 65: GYM STREPTOMYCES MEDIUM

This recipe contains strain-specific modifications for *Nocardiopsis xinjiangensis* DSM 44589 \*

Glucose	4.00	g
Yeast extract	4.00	g
Malt extract	10.00	g
CaCO <sub>3</sub> , if required	2.00	g
Agar	18.00	g
NaCl	150.00	g/l
Distilled water	1000.00	ml

Adjust pH to 7.2 before adding agar. Delete CaCO<sub>3</sub> if liquid medium is used.

\* Plus 15% NaCl