

65: GYM STREPTOMYCES MEDIUM

This recipe contains strain-specific modifications for *Nocardioopsis dassonvillei* DSM 44686 *

| | | |
|---------------------------------|---------|----|
| Glucose | 4.00 | g |
| Yeast extract | 4.00 | g |
| Malt extract | 10.00 | g |
| CaCO ₃ , if required | 2.00 | g |
| Agar | 18.00 | g |
| Distilled water | 1000.00 | ml |

Adjust pH to 7.2 before adding agar. Delete CaCO₃ if liquid medium is used.

* Adjust pH of final medium to 7.2.