



65: GYM STREPTOMYCES MEDIUM

This recipe contains strain-specific modifications for *Nocardiopsis xinjiangensis* DSM 44589 *

Final pH: 7.2

Final volume: 1000 ml

Glucose	4.00	g
Yeast extract	4.00	g
Malt extract	10.00	g
CaCO ₃ , if required	2.00	g
Agar	20.00	g
NaCl	150.00	g/l
Distilled water	1000.00	ml

Adjust pH to 7.2 before adding agar. Delete CaCO₃ if liquid medium is used.

* Plus 15% NaCl