

## 65: GYM STREPTOMYCES MEDIUM

This recipe contains strain-specific modifications for *Saccharomonospora iraqiensis subsp. paurometabolica* DSM 44619 \*

Final pH: 7.2

Final volume: 1000 ml

Glucose	4.00	g
Yeast extract	4.00	g
Malt extract	10.00	g
CaCO <sub>3</sub> , if required	2.00	g
Agar	20.00	g
NaCl	100.00	g/l
Distilled water	1000.00	ml

Adjust pH to 7.2 before adding agar. Delete CaCO<sub>3</sub> if liquid medium is used.

\* Plus 10% NaCl