

65: GYM STREPTOMYCES MEDIUM

This recipe contains strain-specific modifications for *Rugosimonospora acidiphila* DSM 45227 *

Final pH: * 5.5 - 6.0

Final volume: 1000 ml

Glucose	4.00	g
Yeast extract	4.00	g
Malt extract	10.00	g
CaCO ₃ , if required	2.00	g
Agar	20.00	g
Distilled water	1000.00	ml

Adjust pH to 7.2 before adding agar. Delete CaCO₃ if liquid medium is used.

* Adjust pH of final medium to 5.5 - 6.0.