

750: MEDIUM FOR MARINE METHYLOTROPHS

This recipe contains strain-specific modifications for *Methylophaga marina* DSM 5689 *

Final pH: 7.0

Final volume: 1000 ml

Beef extract	7.00	g
Peptone	10.00	g
NaCl	25.00	g
K ₂ HPO ₄	1.00	g
(NH ₄) ₂ SO ₄	1.00	g
Agar	20.00	g
Distilled water	1000.00	ml
Methanol	10.00	ml
MeOH	10.00	g/l

1. pH 7.0
2. After sterilisation, cool to 50°C and add 10 ml methanol.

* Plus 1% MeOH