

830: R2A MEDIUM

This recipe contains strain-specific modifications for *Granulicella arctica* DSM 23128 *

Final pH: * 5.0 - 5.6

Final volume: 1000 ml

Yeast extract	0.50	g
Proteose peptone (Difco no. 3)	0.50	g
Casamino acids	0.50	g
Glucose	0.50	g
Starch (soluble)	0.50	g
Na-pyruvate	0.30	g
K ₂ HPO ₄	0.30	g
MgSO ₄ x 7 H ₂ O	0.05	g
Agar	15.00	g
Agar	5.00	g/l
Distilled water	1000.00	ml

Final pH 7.2; adjust with crystalline K₂HPO₄ or KH₂PO₄ before adding agar. Add agar, heat medium to boiling to dissolve agar, and autoclave for 15 min at 121°C.

* Adjust pH to 5.0 - 5.6 with HCl. Solidify with 5g/l agar