

92: TRYPTICASE SOY YEAST EXTRACT MEDIUM

This recipe contains strain-specific modifications for *Microcella alkalica* DSM 19663 *

Final pH: * 9.0

Final volume: 1000 ml

Trypticase soy broth	30.00	g
Yeast extract	3.00	g
Agar, if necessary	15.00	g
Na-sesquicarbonate	10.00	g/l
Distilled water	1000.00	ml

Adjust pH to 7.0 - 7.2.

* Half-strength plus 1% Na-sesquicarbonate; pH 9.0