

## 92: TRYPTICASE SOY YEAST EXTRACT MEDIUM

This recipe contains strain-specific modifications for *Arthrobacter sp.* DSM 22480 \*

Final pH: \* 7.0 - 8.0

Final volume: 1000 ml

Trypticase soy broth	30.00	g
Yeast extract	3.00	g
Agar, if necessary	15.00	g
Distilled water	1000.00	ml

Adjust pH to 7.0 - 7.2.

\* Adjust pH of final medium to 7.0 - 8.0.