

92: TRYPTICASE SOY YEAST EXTRACT MEDIUM

This recipe contains strain-specific modifications for *Tetragenococcus halophilus subsp. flandriensis* DSM 23766 *

Final pH: 7.0 - 7.2

Final volume: 1000 ml

Trypticase soy broth	30.00	g
Yeast extract	3.00	g
Agar, if necessary	15.00	g
NaCl	50.00	g/l
Distilled water	1000.00	ml

Adjust pH to 7.0 - 7.2.

* Plus NaCl 50 g/L; microaerophilic