

92: TRYPTICASE SOY YEAST EXTRACT MEDIUM

This recipe contains strain-specific modifications for *Intrasporangium sp.* DSM 24459 *

Final pH: 7.0 - 7.2

Final volume: 1000 ml

Trypticase soy broth	30.00	g
Yeast extract	3.00	g
Agar, if necessary	15.00	g
Distilled water	1000.00	ml

Adjust pH to 7.0 - 7.2.

* Plus 100 µl Anaerobiersalze (salt solution of DSMZ Medium 104)