

92: TRYPTICASE SOY YEAST EXTRACT MEDIUM

This recipe contains strain-specific modifications for *Corynebacterium jeikeium* DSM 7112 *

Final pH: 7.0 - 7.2

Final volume: 1000 ml

Trypticase soy broth	30.00	g
Yeast extract	3.00	g
Agar, if necessary	15.00	g
Tween 80	2.00	g/l
Distilled water	1000.00	ml

Adjust pH to 7.0 - 7.2.

* Plus 0.2% tween 80