

993: ISP 5 MEDIUM

This recipe contains strain-specific modifications for *Marinococcus halotolerans* DSM 16375 *

Final pH: 7.0 - 7.4

Final volume: 1000 ml

L-Asparagine (anhydrous basis)	1.00	g
Glycerol	10.00	g
K ₂ HPO ₄ (anhydrous basis)	1.00	g
Distilled water	1000.00	ml
Trace element solution	1.00	ml
Agar	20.00	g
MgCl ₂	100.00	g/l

1. The pH of trace salts solution is about 7.0-7.4. Do not adjust if it is within this range. Liquefy agar by steaming at 100°C for 15-20 minutes.

2. Sterilize in flasks for pouring into Petri dishes. Seven standard 90 mm dishes containing 25 ml per dish.

* Plus 10% MgCl₂

Trace element solution (from medium 84)

FeSO ₄ x 7 H ₂ O	0.10	g
MnCl ₂ x 4 H ₂ O	0.10	g
ZnSO ₄ x 7 H ₂ O	0.10	g
Distilled water	100.00	ml