



1273: 0.01% LB SEAWATER MEDIUM

NaNO ₃	120.00	mg
K ₂ HPO ₄	5.00	mg
Vitamin B ₁₂	1.00	μg
Biotin	1.00	μg
Thiamine HCl	100.00	μg
Fe(III)-EDTA	259.00	μg
Mn(II) x EDTA	332.00	μg
Tris(hydroxymethyl)aminomethane	1.00	g
LB Broth (BD-Difco), Lennox	0.10	g
Sea water	1000.00	ml

1. Mix components thoroughly, adjust pH to 8.0 with HCl and autoclave.
2. Comment: Artificial sea water (such as artificial sea water in Medium No. 1118) can be used instead of natural sea water.