



166: JCM MEDIUM No. 166

Final pH: 7.0

Casamino acids (BD-Difco)	15.00	g
Trisodium citrate	3.00	g
Glutamic acid	2.50	g
MgSO ₄ x 7 H ₂ O	2.50	g
KCl	2.00	g
NaCl	250.00	g
Agar	20.00	g

Add components to distilled water and bring volume to 1.0 L. Adjust pH to 7.0. After autoclaving, readjust pH to 8.5 with a sterile 5% Na₂CO₃ solution.