

Main sol. 634

dehydrated RCM medium (dissolved in distilled water)	1000.00	ml
Sodium resazurin (0.1% w/v)	0.50	ml

Dissolve dehydrated RCM medium (available from BD-BBL, OXOID or MERCK) in distilled water and add 0.5 ml/l Na-resazurin solution (0.1% w/v), bring to the boil and cool to room temperature while sparging with 100% N₂ gas. Dispense medium under same gas atmosphere into anoxic Hungate-type tubes or serum vials and autoclave.