

Solution A

Brain heart infusion broth	12.50	g
Trypticase (BBL 11921)	10.00	g
Yeast extract	2.50	g
Na-thioglycolate	0.50	g
L-Cysteine HCl	1.00	g
L-Asparagine	0.25	g
Glucose	2.00	g
Agar	3.00	g
Distilled water	1000.00	ml

Boil for 5 min, cool to 50°C while gassing with 80% N₂ + 20% CO₂. Distribute to anaerobic tubes under the same atmosphere and autoclave.