



Main sol. 121

Na-aspartate	10.00	g
Yeast extract	0.20	g
MgSO ₄ x 7 H ₂ O	1.00	g
CaCl ₂ x 2 H ₂ O	28.00	mg
K ₂ HPO ₄	0.75	g
NaH ₂ PO ₄	0.25	g
L-Cysteine HCl x H ₂ O	0.25	g
Resazurin	1.00	mg
Trace element solution SL-10	1.00	ml
Distilled water	1000.00	ml

Adjust pH to 7.0. Phosphate and cysteine should be dissolved together in 1/10 vol. of water and sterilized separately.