

Main sol. 125 (from medium 125)

KNO ₃	1.00	g
MgSO ₄ x 7 H ₂ O	0.20	g
CaCl ₂ x 2 H ₂ O	0.02	g
Na ₂ HPO ₄	0.23	g
NaH ₂ PO ₄	0.07	g
FeSO ₄ x 7 H ₂ O	1.00	mg
CuSO ₄ x 5 H ₂ O	5.00	µg
H ₃ BO ₃	10.00	µg
MnSO ₄ x 5 H ₂ O	10.00	µg
ZnSO ₄ x 7 H ₂ O	70.00	µg
MoO ₃	10.00	µg
Agar	12.00	g
Distilled water	1000.00	ml

Add 5 ml methanol or incubate under methane-air (4:1). Adjust pH to 6.8.