



SSE (double concentrated)

| | | |
|--|-----------|----|
| CaCl ₂ x 2 H ₂ O | 0.2938 | g |
| NH ₄ Cl | 0.1069 | g |
| MgCl ₂ x 6 H ₂ O | 0.2036 | g |
| (NH ₄) ₂ SO ₄ | 0.1983 | g |
| MgSO ₄ x 7 H ₂ O | 0.7390 | g |
| CaSO ₄ x 2 H ₂ O | 0.8606 | g |
| Ca(NO ₃) ₂ x 4 H ₂ O | 0.2360 | g |
| Na NO ₃ | 0.4240 | g |
| KH ₂ PO ₄ | 0.0700 | g |
| FeSO ₄ x 7 H ₂ O | 0.0111 | g |
| K ₂ SO ₄ | 0.0870 | g |
| Distilled water | 1000.0000 | ml |

Dissolve by shaking overnight. Sterilize by autoclaving.