

## Solution C

Calcium D-(+)-pantothenate	10.00	mg
Choline chloride	100.00	mg
Myo-inositol	100.00	mg
Thiamine HCl	10.00	mg
Thiamin pyrophosphate	10.00	mg
Coenzym A, sodium salt	2.00	mg
Flavin adenine dinucleotide (FAD)	2.00	mg
β-nicotinamide adenine dinucleotide (NAD)	2.00	mg
2-Mercaptoethanesulfonic acid (MES), sodium salt	20.00	mg
N-acetylmuramic acid	50.00	mg
N-Acetylglucosamine	400.00	mg
D-glucuronic acid, sodium salt	1600.00	mg
D-galacturonic acid	1600.00	mg
Pyruvic acid, sodium salt	1100.00	mg
Fumaric acid, disodium salt	1000.00	mg
Formic acid, sodium salt	600.00	mg
Distilled water	100.00	ml

Prepare Solutions C, flush with nitrogen to make them anoxic, filter-sterilize and store at -20°C