

10: ZYMOMONAS MEDIUM

This recipe contains strain-specific modifications for *Zymomonas mobilis* DSM 424 *

Final volume: 1000 ml

Bacto peptone	10.00	g
Yeast extract	10.00	g
Glucose	20.00	g
Agar, for solid medium	15.00	g
Distilled water	1000.00	ml

Incubate microaerobic, for example in a candle jar.

* Anaerobic. Rehydrate and grow lyophilized cells from the ampoule in liquid medium 10. Subsequent subculturing may be carried out in liquid medium or on agar plates.