## **Microorganisms**



## 1: NUTRIENT AGAR

This recipe contains strain-specific modifications for Oxalicibacterium flavum DSM 15506 \*

Final pH: 7.0

Final volume: 1000 ml

Peptone	5.00	g
Meat extract	3.00	g
Agar, for solid medium	15.00	g
Glycerol	10.00	g/l
Na-oxalate	4.00	g/l
Distilled water	1000.00	ml

- 1. Adjust pH to 7.0.
- 2. For Bacillus strains the addition of 10.0 mg MnSO $_4$  x  $H_2O$  is recommended for sporulation.

<sup>\*</sup> With 1% glycerol and 0.4% K- or Na-oxalate; rehydrate and grow lyophylisate in liquid medium 1. Subsequently, subculturing may be carried out in liquid medium or on agar plates