

## 535: TRYPTICASE SOY BROTH AGAR

This recipe contains strain-specific modifications for *Streptobacillus ratti* DSM 101843 \*

Final pH: 7.3

Final volume: 1000 ml

|                                                              |         |     |
|--------------------------------------------------------------|---------|-----|
| Trypticase soy broth (BBL 11768, Oxoid CM129 or Merck 30460) | 30.00   | g   |
| Agar                                                         | 20.00   | g   |
| Horse serum                                                  | 200.00  | g/l |
| Distilled water                                              | 1000.00 | ml  |

1. pH 7.3
2. Autoclave at 121°C for 15 min.

\* Plus 20% horse serum; microaerophilic