

535: TRYPTICASE SOY BROTH AGAR

This recipe contains strain-specific modifications for *Nocardioides daejeonensis* DSM 25320 *

Final pH: 7.3

Final volume: 1000 ml

Trypticase soy broth (BBL 11768, Oxoid CM129 or Merck 30460)	30.00	g
Agar	20.00	g
Tween 80	1.00	g/l
Distilled water	1000.00	ml

1. pH 7.3
2. Autoclave at 121°C for 15 min.

* Plus 0,1% Tween 80