

535: TRYPTICASE SOY BROTH AGAR

This recipe contains strain-specific modifications for *Dermatophilus congolensis* DSM 44174 *

Final pH: 7.3

Final volume: 1000 ml

Trypticase soy broth (BBL 11768, Oxoid CM129 or Merck 30460)	30.00	g
Agar	20.00	g
Blood	5.00	%(v/v)
Distilled water	1000.00	ml

1. pH 7.3
2. Autoclave at 121°C for 15 min.

* With 5% (v/v) blood