

## 535: TRYPTICASE SOY BROTH AGAR

This recipe contains strain-specific modifications for *Corynebacterium coyleae* DSM 44185 \*

Final pH: 7.3

Final volume: 1000 ml

|                                                              |         |        |
|--------------------------------------------------------------|---------|--------|
| Trypticase soy broth (BBL 11768, Oxoid CM129 or Merck 30460) | 30.00   | g      |
| Agar                                                         | 20.00   | g      |
| Blood                                                        | 5.00    | %(v/v) |
| Distilled water                                              | 1000.00 | ml     |

1. pH 7.3
2. Autoclave at 121°C for 15 min.

\* By choice with 5% (v/v) blood; obligately anaerobic