

## 535: TRYPTICASE SOY BROTH AGAR

This recipe contains strain-specific modifications for *Cephalotococcus capnophilus* DSM 100879 \*

Final pH: 7.3

Final volume: 1000 ml

Trypticase soy broth (BBL 11768, Oxoid CM129 or Merck 30469)	30.00	g
Agar	20.00	g
Distilled water	1000.00	ml

1. pH 7.3
2. Autoclave at 121°C for 15 min.

\* Microaerophilic under 93% N<sub>2</sub>, 5% CO<sub>2</sub>, 2% O<sub>2</sub>