

535: TRYPTICASE SOY BROTH AGAR

This recipe contains strain-specific modifications for *Nocardia asteroides* DSM 44145 *

Final pH: 7.3

Final volume: 1000 ml

| | | |
|--|---------|-----|
| Trypticase soy broth (BBL 11768, Oxoid CM129 or Merck 30460) | 30.00 | g |
| Agar | 20.00 | g |
| Defibrinated sheep blood | 50.00 | g/l |
| Distilled water | 1000.00 | ml |

1. pH 7.3
2. Autoclave at 121°C for 15 min.

* With 5% defibrinated sheep blood