



## 535: TRYPTICASE SOY BROTH AGAR

This recipe contains strain-specific modifications for *Corynebacterium coyleae* DSM 44185 \*

Final pH: 7.3

Final volume: 1000 ml

Trypticase soy broth (BBL 11768, Oxoid CM129 or Merck <del>30469</del> )	30.00	g
Agar	20.00	g
Blood	5.00	%(v/v)
Distilled water	1000.00	ml

1. pH 7.3
2. Autoclave at 121°C for 15 min.

\* By choice with 5% (v/v) blood; obligately anaerobic