

661: DP MEDIUM

This recipe contains strain-specific modifications for *Escherichia coli* DSM 8590 *

Final volume: 1000 ml

Tryptone	10.00	g
Yeast extract	10.00	g
Glucose	5.00	g
KH ₂ PO ₄	3.00	g
K ₂ HPO ₄	7.00	g
Na ₃ -citrate x 3 H ₂ O	0.50	g
MgSO ₄ x 7 H ₂ O	0.10	g
Agar	2.00	%(w/v)
Erythromycin	50.00	µg/ml
Distilled water	1000.00	ml

* Plus 50 µg/ml erythromycin; Please, see hints on growth of this strain