

827: THERMOANAEROBACTER SULFUROPHILUS MEDIUM

This recipe contains strain-specific modifications for *Petrimonas sulfuriphila* DSM 16547 *

Final pH: 6.8 - 7.0

Final volume: 1002 ml

NH ₄ Cl	0.33	g
KCl	0.33	g
KH ₂ PO ₄	0.33	g
MgCl ₂ x 6 H ₂ O	0.33	g
CaCl ₂ x 2 H ₂ O	0.33	g
Trace element solution SL-10	1.00	ml
Sodium resazurin (0.1% w/v)	0.50	ml
Sulfur (powdered)	10.00	g
Na ₂ CO ₃	1.00	g
Yeast extract	0.50	g
D-Glucose	5.00	g
Wolin's vitamin solution (10x)	1.00	ml
Na ₂ S x 9 H ₂ O	0.50	g
Trypticase peptone	3.00	g
Distilled water	1000.00	ml

Dissolve ingredients except sulfur, carbonate, yeast extract, glucose, vitamins, and sulfide, then sparge medium with 80% N₂ and 20% CO₂ gas mixture for 30 - 45 min to make it anoxic. Dispense medium under same gas atmosphere into anoxic Hungate-type tubes or serum vials containing already the appropriate amount of sulfur. Heat the medium vials to 90-100°C for 1 hour on each of 3 successive days to sterilize the medium. Add yeast extract, glucose, vitamins (sterilized by filtration), and sulfide from sterile anoxic stock solutions prepared under 100% N₂ gas and carbonate from a sterile anoxic stock solution prepared under 80% N₂ and 20% CO₂ gas mixture. Prior to use check pH of complete medium and adjust to 6.8 - 7.0, if necessary.

* Supplement medium with 3.00 g/l Trypticase peptone.

Trace element solution SL-10 (from medium 320)

HCl (25%)	10.00	ml
FeCl ₂ x 4 H ₂ O	1.50	g
ZnCl ₂	70.00	mg
MnCl ₂ x 4 H ₂ O	100.00	mg
H ₃ BO ₃	6.00	mg
CoCl ₂ x 6 H ₂ O	190.00	mg
CuCl ₂ x 2 H ₂ O	2.00	mg
NiCl ₂ x 6 H ₂ O	24.00	mg

827: THERMOANAEROBACTER SULFUROPHILUS MEDIUM

Na ₂ MoO ₄ x 2 H ₂ O	36.00	mg
Distilled water	990.00	ml

First dissolve FeCl₂ in the HCl, then dilute in water, add and dissolve the other salts. Finally make up to 1000.00 ml.

Wolin's vitamin solution (10x) (from medium 120)

Biotin	20.00	mg
Folic acid	20.00	mg
Pyridoxine hydrochloride	100.00	mg
Thiamine HCl	50.00	mg
Riboflavin	50.00	mg
Nicotinic acid	50.00	mg
Calcium D-(+)-pantothenate	50.00	mg
Vitamin B ₁₂	1.00	mg
p-Aminobenzoic acid	50.00	mg
(DL)-alpha-Lipoic acid	50.00	mg
Distilled water	1000.00	ml