Microorganisms



838: GEOANAEROBACTER MEDIUM

This recipe contains strain-specific modifications for $Geopsychrobacter\ electrodiphilus\ DSM\ 16401\ ^*$

Final pH: 6.8

Final volume: 1002 ml

Ferric citrate monohydrate	10.00	g	
KH ₂ PO ₄	0.60	g	
NH ₄ Cl	0.30	g	
$MgSO_4 \times 7 H_2O$	0.50	g	
CaCl ₂ x 2 H ₂ O	0.10	g	
Trace element solution SL-11	1.00	ml	
Selenite-tungstate solution	1.00	ml	
Na ₂ CO ₃	1.50	g	
Na-acetate	0.80	g	
Na ascorbate	0.80	g	
NaCl	20.00	_	
IVaCI	20.00	g	
Wolin's vitamin solution	10.00	ml	
Distilled water	1000.00	ml	

First dissolve ferric citrate by heating the water under continuous stirring. After cooling to room temperature adjust the pH to 6.0, then add and dissolve the remaining ingredients, except carbonate, acetate and ascorbate. Sparge medium with 80% N_2 and 20% CO_2 gas mixture for 30 - 45 min to make it anoxic, then dispense under the same gas atmosphere into anoxic Hungate-type tubes or serum vials and autoclave. After autoclaving add acetate and ascorbate from sterile anoxic stock solutions prepared under 100% N_2 gas and carbonate from a sterile anoxic stock solution prepared under 80% N_2 and 20% CO_2 gas mixture. Prior to inoculation check the medium pH and adjust to 6.8, if necessary.

Trace element solution SL-11 (from medium 722)

	•	
Na_2 -EDTA x 2 H_2O	5.20	g
FeCl ₂ x 4 H ₂ O	1.50	g
ZnCl ₂	70.00	mg
$MnCl_2 \times 4 H_2O$	100.00	mg
H ₃ BO ₃	6.00	mg
CoCl ₂ x 6 H ₂ O	190.00	mg
CuCl ₂ x 2 H ₂ O	2.00	mg
$NiCl_2 \times 6 H_2O$	24.00	mg

^{*} Supplement medium with 20.00 g/l NaCl and omit Na-ascorbate. After autoclaving supplement medium with 10.00 ml/l of a sterile anoxic Wolin's vitamin solution (see medium 141).

Microorganisms

838: GEOANAEROBACTER MEDIUM



 $Na_2MoO_4 \times 2 H_2O$ 36.00 mg Distilled water 1000.00 ml

Dissolve EDTA in 800 ml distilled water, adjust pH to 7 using 2 N NaOH and add ferrous chloride. After ferrous chloride has dissolved add other compounds. Finally adjust pH to 6.0 and bring volume to 1000 ml.

Selenite-tungstate solution (from medium 385)

NaOH	0.50	g
$Na_2SeO_3 \times 5 H_2O$	3.00	mg
$Na_2WO_4 \times 2 H_2O$	4.00	mg
Distilled water	1000.00	ml

Wolin's vitamin solution* (from medium 141)

Biotin	2.00	mg
Folic acid	2.00	mg
Pyridoxine hydrochloride	10.00	mg
Thiamine HCl	5.00	mg
Riboflavin	5.00	mg
Nicotinic acid	5.00	mg
Calcium D-(+)-pantothenate	5.00	mg
Vitamin B ₁₂	0.10	mg
p-Aminobenzoic acid	5.00	mg
(DL)-alpha-Lipoic acid	5.00	mg
Distilled water	1000.00	ml